

Clinical and Biophysical Efficacy of Olav Skille's Original Vibroacoustic Therapy in Osteoporosis and Bone Density Rehabilitation

Historical Genesis and Evolution of Olav Skille's Vibroacoustic Therapy

Vibroacoustic Therapy (VAT) was developed in Norway in the late 1960s and early 1980s by the psychologist, therapist, and musician Olav Skille.¹ Skille's initial exploration began in 1968 with the creation of the Musical Behavior Scale (MUBS), a nonverbal screening instrument designed to evaluate sensory, motor, and psychological parameters in children with severe developmental and cognitive disabilities.³ Through close observation of how physical sound waves interacted with these patients, Skille and his collaborator Juliette Alvin recognized that audible music could elicit profound physiological changes that surpassed mere emotional or psychological responses.³ This realization led Skille to isolate pure, low-frequency sinusoidal sound waves and apply them directly to the human body, officially launching the clinical discipline of Vibroacoustic Therapy in 1980.⁴

By 1982, Skille established the precise technical definition of VAT as the therapeutic application of low-frequency sinusoidal sound waves—constrained within the range of **30 Hz** to **120 Hz**—delivered in tandem with specially composed, rhythmically neutral music.³ In 1987, Skille published his first seminal paper documenting the multi-center clinical efficacy of VAT across Germany, Finland, and England, demonstrating its utility in treating muscle spasms, spasticity in cerebral palsy, asthma, and chronic pain.⁴

A critical distinction must be maintained between original Vibroacoustic Therapy and "vibroacoustic music" (VAM).⁴ VAM typically involves streaming commercial music overlaid with a generic low-frequency hum.⁴ Conversely, the original Skille method, which forms the bedrock of modern clinical VAT, relies strictly on single, mathematically precise sinusoidal waves.⁴ Pure sine waves contain no overtones, harmonic distortions, or acoustic transients, allowing them to target specific tissues without inducing sensory fatigue or neural noise.⁴

Throughout his career, which spanned over 30,000 hours of documented clinical studies, Skille researched the specific therapeutic windows of different frequencies.⁸ This legacy of clinical inquiry is preserved and advanced in North America by TheSoundWell, under the leadership of CEO Avigail Berg, who was personally trained by Olav Skille.⁴ TheSoundWell systems utilize Skille's original therapeutic frequencies, integrating his clinical methods with modern tactile

sound technology to address chronic musculoskeletal, vascular, and neurological disorders.¹⁰

Biophysical Foundations of Low-Frequency Sound and Human Tissue Interaction

The therapeutic efficacy of Vibroacoustic Therapy is rooted in the physics of acoustic propagation through liquid media.⁴ Because the human body is composed of approximately 60% to 70% water, it acts as an exceptional conductor of acoustic waves.⁴ Acoustic waves travel approximately five times faster and with significantly less attenuation through human soft tissue and bone than through air.⁴ This tissue conductivity is governed by the classical wave equation:

$$v = f \cdot \lambda$$

where v is the velocity of sound in human tissue (approximately 1500 m/s), f is the frequency of the wave, and λ is the wavelength.⁴ A 40 Hz pure sine wave applied directly to the body produces a wavelength of approximately 37.5 meters .⁴

This long wavelength allows the acoustic energy to easily bypass superficial epidermal and fascial resistance, penetrating deep into the visceral organs, nervous system, and skeletal structures.⁴ This deep penetration creates a phenomenon known as sympathetic resonance, wherein the physical tissues of the body vibrate in harmony with the external acoustic stimulus, acting as a deep cellular-level massage.⁴

The physical transfer of this acoustic energy from digital files to the human skeleton requires specialized hardware.¹ TheSoundWell equipment, including the Vibroacoustic Therapy Vibrating Mat, utilizes six electromagnetic transducers embedded within high-density, acoustically receptive foam.⁴ These transducers act as specialized speakers that convert the electrical output of a frequency generator into physical, non-turbulent sinusoidal oscillations.¹

When a patient lies in a supine position on the mat, this low-frequency sound stimulation (LFSS) initiates a predictable, time-dependent physiological cascade:

- **Endothelial Stimulation (Within 1 Minute):** The mechanical shear stress of the acoustic wave stimulates the endothelial cells lining the microvasculature, triggering an immediate release of nitric oxide (NO) and Adrenomedullin (AM).⁴ Nitric oxide relaxes the smooth muscle of blood vessel walls, accelerating local blood flow, enhancing oxygenation, and reducing deep-tissue inflammation.⁴
- **Mechanoreceptor Activation (Within 3 Minutes):** The physical oscillations stimulate low-threshold mechanoreceptors within the skin and deep fascial tissues, particularly the Pacinian corpuscles.¹ These receptors send afferent signals through the vagus nerve to

the brainstem, shifting the autonomic nervous system from sympathetic dominance to parasympathetic regulation.⁴

- **Systemic Homeostasis (Within 10 Minutes):** Continuous harmonic stimulation drives the nervous system into a balanced state.¹⁰ Cortisol levels decrease, while the production of serotonin and endorphins increases, allowing muscle tension to dissolve and supporting systemic tissue recovery.¹⁰

Mechanotransduction and Cellular Pathways of Bone Remodeling

Bone tissue is a highly dynamic organ that constantly adapts its micro-architecture to external mechanical demands.¹⁸ This adaptation is managed within bone multicellular units (BMUs) through bone remodeling—a continuous process of bone resorption by osteoclasts followed by bone deposition and mineralization by osteoblasts.²⁰ In patients with osteopenia and osteoporosis, this balance is disrupted; osteoclastic activity outpaces osteoblastic bone formation, resulting in the thinning of trabecular struts, increased skeletal fragility, and an elevated risk of fracture.²¹

Low-frequency sound and mechanical vibration serve as potent physical stimuli that can restore bone homeostasis by initiating cellular mechanotransduction.¹⁸


At the cellular scale, osteocytes—which reside within the lacunar-canalicular network—act as the primary mechanosensors of the skeleton.¹⁸ When low-frequency acoustic vibrations pass through bone, they generate microscopic pressure waves that cause interstitial fluid to flow through these canaliculi.¹⁸ This fluid flow exerts Newtonian fluid shear stress on the cell membranes of both osteocytes and bone marrow stromal cells (BMSCs).¹⁸

The physical shear stress (τ) is mathematically defined as:

$$\tau = \mu \cdot \frac{du}{dy}$$

where μ represents the dynamic viscosity of the interstitial fluid and $\frac{du}{dy}$ represents the fluid velocity gradient across the cell surface.²³ This mechanical force is detected by primary cilia, hair-like sensory organelles projecting from the bone cell membrane.²³

Deflection of these primary cilia opens stretch-activated calcium channels, inducing an influx of extracellular calcium (Ca^{2+}) into the cytoplasm.²³ This intracellular calcium signal triggers a signaling cascade that activates the *Wnt*/ β -catenin pathway.¹⁸



The activation of the *Wnt/β-catenin* pathway upregulates the transcription factors Runt-related transcription factor 2 (Runx2) and Osterix (OSX).¹⁸ Runx2 and Osterix direct the differentiation of multipotent BMSCs down the osteogenic lineage, promoting their maturation into active, osteoid-secreting osteoblasts.²³ These mature osteoblasts synthesize type I collagen and secrete bone morphogenetic protein 2 (BMP-2), accelerating matrix deposition and subsequent mineralization.¹⁸

Simultaneously, mechanical stimulation upregulates the expression of osteoprotegerin (OPG)

and downregulates receptor activator of nuclear factor κ B ligand (RANKL) in osteoblastic cells.¹⁸ By decreasing the RANKL/OPG ratio, OPG competitively binds to RANK on osteoclast

precursors, blocking their fusion and activation, which halts osteoclast-mediated bone demineralization.¹⁸

Skeletal tissues exhibit a distinct biophysical resonance window between **25 Hz** and **50 Hz**.¹⁰ In vitro and in vivo studies indicate that stimulation within this range—specifically at **30 Hz**, **45 Hz**, and **50 Hz**—elicits the most pronounced osteogenic gene expression and osteoblast proliferation.²³ This biological resonance corresponds to the lower frequency range of Olav Skille's Vibroacoustic Therapy, which explicitly designates the **25–50 Hz** range for bone healing and density support.¹⁰

Clinical Evidence of Vibration and Sound-Vibration Efficacy in Osteoporosis

While clinical research directly evaluating Olav Skille's acoustic platforms for osteoporosis is ongoing, the clinical efficacy of low-frequency, low-magnitude physical vibrations of identical frequency profiles is well-documented.¹⁴ A pivotal study by Dr. Lee Bartel and Abdullah Mosabbir at the University of Toronto demonstrated that the human body does not differentiate between mechanical and acoustic vibrations at these low frequencies.¹⁶

Both modalities activate the same physiological, hemodynamic, and skeletal pathways.¹⁶ Consequently, clinical trials utilizing low-intensity mechanical vibrations of identical frequencies (**30–50 Hz**) provide strong support for the clinical utility of acoustic vibroacoustic therapy.¹⁴

In a randomized, double-blind, placebo-controlled trial of 70 postmenopausal women, Rubin et al. evaluated the effects of daily low-intensity vibration (LIV) applied at **30 Hz** with an acceleration of **0.2g** (**2 m/s²**) for 20 minutes daily over 12 months.²³ The active vibration group experienced a 1.5% reduction in bone loss in the spine and a 2.17% reduction in the femoral neck compared to the control group.²⁶ Notably, compliant women weighing under **65 kg** in the active vibration group achieved an average increase in bone mineral density (BMD) of 3.35%.²⁶

Similarly, Ruan et al. conducted a 6-month randomized controlled trial involving 116 postmenopausal women with established osteoporosis.²⁶ The treatment group received **30 Hz** vibrations with a **5 mm** amplitude for 10 minutes, five times per week.²⁶ At the end of 6 months, the control group showed a 1.9% decrease in femoral neck BMD, while the vibration group achieved a 3.2% increase in femoral neck BMD, demonstrating that low-frequency vibration can help reverse osteoporotic bone loss.²⁶

Further clinical support for targeted vibration comes from the FDA-cleared Osteoboost belt, a

wearable device designed to prevent osteopenia and osteoporosis using low-magnitude high-frequency vibration (LMHFV).²⁷ In a 12-month randomized, double-blind, sham-controlled trial of 126 postmenopausal women, participants used the active or sham device for 30 minutes per day, at least five times per week.²⁸

The active group demonstrated an 85% reduction in vertebral bone density loss, an 83% reduction in vertebral bone strength loss, and a 55% reduction in hip bone density loss compared to the sham group, with zero device-related serious adverse events.²¹

These findings provide strong clinical evidence that targeted, low-intensity vibration at appropriate frequencies can safely and effectively preserve bone density and strength.²⁷

Study / Source	Clinical Model & Population	Device Type & Parameter Profile	Key Quantitative & Clinical Outcomes
Rubin et al. ²³	70 postmenopausal women; 12-month double-blind, RCT.	Standing low-intensity platform; 30 Hz frequency, 0.2g acceleration, 2 × . ²³	Reduced spinal bone loss by 1.5% and femoral neck loss by 2.17%; compliant women < gained 3.35% BMD. ²⁶
Ruan et al. ²⁶	116 postmenopausal women with osteoporosis; 6-month RCT.	Standing vibration platform; 30 Hz frequency, 5 mm amplitude, 10 min/day, 5 days/week. ²⁶	Control group BMD decreased by 1.9%; treatment group femoral neck BMD increased by 3.2%. ²⁶
Osteoboost Trial ²⁷	126 postmenopausal women with low bone density; 12-month sham-controlled RCT.	Wearable vibration belt; targeted mechanical stimulation to hips/spine, 30 min, 5 days/week. ²⁷	85% reduction in vertebral bone density loss, 83% reduction in vertebral strength loss, and 55% reduction in hip bone density loss. ²¹

<p>BMSC In Vivo Trial ²⁵</p>	<p>New Zealand rabbits with bone injuries; comparative animal study.</p>	<p>Experimental vibration platform; compared 12.5, 25, 50, 100, and 200 Hz frequencies.²⁵</p>	<p>Identified 50 Hz as the optimal frequency to stimulate BMSC differentiation, showing faster callus absorption and orderly collagen alignment.²⁵</p>
<p>Broad Frequency Trial¹⁹</p>	<p>In vivo mouse ulnar loading model; comparative experimental trial.</p>	<p>Combined ulnar loading (2 Hz) with broad-frequency vibration (0–50 Hz), 3 min/day.¹⁹</p>	<p>Combined loading enhanced periosteal bone formation by 3.9-fold compared to standard exercise loading alone.¹⁹</p>

Strategic Application of TheSoundWell Ecosystem in Osteoporosis Management

Integrating Olav Skille’s Vibroacoustic Therapy into the comprehensive care of osteoporosis patients offers several clinical advantages over conventional physical interventions.⁴ Patients with severe osteopenia or osteoporosis are often unable to engage in high-impact weight-bearing exercises due to pain, physical weakness, or an increased risk of fall-induced fractures.²²


TheSoundWell equipment, including the Vibroacoustic Mat and Recliner, allows patients to receive osteogenic stimulation in a safe, fully supported, non-weight-bearing supine position.⁴ This setup eliminates gravitational load on fragile joints and fragile vertebrae, providing an accessible and comfortable therapy option.⁴

Furthermore, the biophysical concept of stochastic resonance can be utilized to enhance bone remodeling.¹⁹ Stochastic resonance occurs when a weak, sub-threshold sensory signal is amplified by the addition of low-amplitude broad-frequency noise.¹⁹

In vivo bone research demonstrates that while passive, low-amplitude vibration alone does not trigger new bone formation, combining a broad-frequency vibration (**0–50 Hz**) with standard physical exercise enhances the osteogenic response on the bone surface by nearly 4-fold (3.9x) compared to exercise alone.¹⁹

By administering TheSoundWell's low-amplitude sinusoidal vibrations, clinicians can lower the

mechanical threshold required to trigger bone remodeling, significantly boosting the bone-building efficacy of light physical exercise or standard physiotherapy.¹⁰



Additionally, VAT helps mitigate stress-induced bone loss.¹⁰ Chronic anxiety, PTSD, and severe psychological distress trigger hyperactivation of the hypothalamic-pituitary-adrenal (HPA) axis, leading to chronically elevated levels of systemic cortisol.¹⁰ High cortisol levels suppress osteoblast differentiation and survival while upregulating RANKL, accelerating glucocorticoid-induced bone loss.¹⁰

Because TheSoundWell VAT systems are designed to shift the autonomic nervous system to parasympathetic dominance and reduce cortisol levels, they help address this endocrine pathway of bone loss.¹⁰ By lowering cortisol and restoring HPA axis homeostasis, VAT helps shift the systemic environment from a catabolic state to an anabolic state, supporting natural bone preservation.¹⁰

Parameter	Vibroacoustic Therapy (VAT)	Whole-Body Vibration (WBV) [Mechanical Platforms]
Primary Physical Stimulus	Pure sinusoidal sound waves via electro-magnetic transducers ¹	Motor-driven mechanical oscillations via eccentric flywheels ²²
Frequency Profile	Precise 30–120 H sine waves; zero mechanical harmonics ³	5–90 H oscillations with variable mechanical noise ⁷
Mechanical Intensity	Micro-displacement; extremely low acceleration magnitude ($< 0.001 g$) ⁴	Macro-displacement (0.5–10 mm); high acceleration (0.3 g to $> 1 g$) ⁷
Patient Positioning	Supine or reclined; fully supported on ergonomic mats/recliners ⁴	Upright standing (bipedal/unipedal) or active semi-squat positions ²²
Biomechanical Loading	Non-weight-bearing; passive tissue and fluid resonance ⁴	Active muscle contractions; compressive joint and spinal loading ²²
Skeletal Risk Profile	Safe; zero joint friction or impact-induced spinal micro-fracture risks ¹⁰	Moderate-high; standing impact may cause pain or damage in fragile bones ¹⁸
Fall & Balance Hazard	Non-existent; ideal for patients with advanced ataxia or weakness ⁴	Present; high-intensity oscillation can induce dizziness or balance loss ²²
Endocrine Response	Suppresses HPA axis, lowers cortisol, increases endorphins ¹⁰	Stimulates HPA axis transiently; increases serum GH/testosterone ⁷
Vasodilatory Mechanism	Rapid endothelial nitric oxide (NO) release ⁴	Muscle contraction-driven mechanical pumping ²²

User Compliance	High; comfortable, passive, soothing deep-tissue massage ⁴	Moderate; physically taxing, can cause joint discomfort in some users ²²
------------------------	---	---

Clinical Protocols, Contraindications, and Safety Standards

To safely integrate Vibroacoustic Therapy into an osteoporosis rehabilitation program, clinicians should follow structured protocols and establish clear safety standards.¹⁰ TheSoundWell equipment should be configured to target skeletal structures while ensuring patient comfort and compliance ⁴:

- Acoustic Frequency Parameters:** Clinicians should select Skille's foundational sinusoidal frequencies within the **25–50 Hz** range, which targets bone tissue resonance and BMSC mechanotransduction.¹⁰ Frequencies above this range, such as **40 Hz**, can be included to promote relaxation and neural coherence, while lower frequencies directly address local bone remodeling.¹
- Session Duration and Frequency:** A standard therapeutic session should run for 23 to 30 minutes, 3 to 5 times per week.¹⁰ The effects of low-frequency sound on bone tissue are cumulative, requiring consistent application over a minimum of 6 to 12 months to achieve measurable changes in bone mineral density.²⁶
- Equipment and Comfort Setup:** The patient is placed on TheSoundWell Vibroacoustic Mat or Recliner in a supine position.⁴ For patients with severe spinal kyphosis or chronic back pain, conforming support systems such as the Yogibo Max Body Pillow should be used to minimize pressure points and ensure uniform acoustic wave transmission.¹⁵
- Systemic Interventions:** To support bone mineralization, patients should maintain standardized daily intake of calcium and vitamin D.²⁸ Combining VAT with light, safe resistance or balance exercises can also leverage stochastic resonance to optimize the osteogenic response.¹⁹

While Vibroacoustic Therapy is a safe and gentle modality, clinicians must observe specific precautions and contraindications ¹⁰:

- Absolute Contraindications:** VAT must not be used in patients with active deep vein thrombosis (DVT) or thrombophlebitis, as acoustic-induced vasodilation and accelerated blood flow could increase thromboembolic risks.¹⁰ It is also contraindicated for patients with unhealed, unstable skeletal fractures, severe epilepsy triggered by sound or vibration, or implanted electronic devices (such as pacemakers or defibrillators) unless cleared by their cardiologist.¹⁰
- Relative Precautions:** Patients with severe, fragile osteoporosis or active bone tumors should consult their primary physician before starting treatment.¹⁰ During sessions, the

physical vibration intensity must always be kept within the patient's comfortable range, avoiding any harsh or jarring sensations.⁴ Proper hydration before and after sessions is also recommended to support optimal vascular and tissue resonance.³⁶

By adhering to these clinical standards, healthcare providers can safely integrate TheSoundWell's vibroacoustic systems into multidisciplinary osteoporosis and osteopenia care, offering a gentle, non-weight-bearing option that supports bone preservation and overall patient well-being.⁴

Works cited

1. Vibroacoustic therapy - Wikipedia, accessed May 19, 2026, https://en.wikipedia.org/wiki/Vibroacoustic_therapy
2. Vibroacoustic Therapy: Guide to Healing w/ Sound & Vibration - Carbon Wellness, accessed May 19, 2026, <https://carbonwellnessmd.com/blogs/vibroacoustic/vibroacoustic-therapy-healing>
3. VibroAcoustic Therapy - Oxford Academic, accessed May 19, 2026, <https://academic.oup.com/musictherapy/article-pdf/8/1/61/8739474/8-1-61.pdf>
4. The Magic of Original Vibroacoustic Therapy: How Olav Skille's Method, Brought to Life by TheSoundWell, Helps Neurodivergent Kids and Their Families Chill Out and Thrive | by Avigaili Berg | Medium, accessed May 19, 2026, <https://medium.com/@gaili/the-magic-of-original-vibroacoustic-therapy-how-olav-s-killies-method-brought-to-life-by-66a272bf731d>
5. Therapeutic Modalities - Non-Thermal | PM&R KnowledgeNow - AAPM&R, accessed May 19, 2026, <https://now.aapmr.org/therapeutic-modalities/>
6. Effects of Vibroacoustic Stimulation on Psychological, Physiological, and Cognitive Stress, accessed May 19, 2026, <https://pmc.ncbi.nlm.nih.gov/articles/PMC11436230/>
7. Potential of Vibroacoustic Therapy in Persons with Cerebral Palsy: An Advanced Narrative Review - PMC, accessed May 19, 2026, <https://pmc.ncbi.nlm.nih.gov/articles/PMC6843730/>
8. vibroacoustic therapy: conception, development, and future directions - ResearchGate, accessed May 19, 2026, https://www.researchgate.net/publication/332947314_VIBROACOUSTIC_THERAPY_CONCEPTION_DEVELOPMENT_AND_FUTURE_DIRECTIONS
9. VIBROACOUSTIC THERAPY: CONCEPTION, DEVELOPMENT, AND FUTURE DIRECTIONS - ResearchGate, accessed May 19, 2026, https://www.researchgate.net/profile/Elsa-Campbell/publication/332947314_VIBROACOUSTIC_THERAPY_CONCEPTION_DEVELOPMENT_AND_FUTURE_DIRECTIONS/links/5cd2e81c299bf14d9581698a/VIBROACOUSTIC-THERAPY-CONCEPTION-DEVELOPMENT-AND-FUTURE-DIRECTIONS.pdf
10. TheSoundWell Vibro-therapy | anxiety support network | 3035 ..., accessed May 19, 2026, <https://www.vibro-therapy.com/>
11. Reiki Sound Healing Clermont | Music Therapy Classes, Sound Bath, accessed May 19, 2026, <https://aunaturalespa.com/sound-healing/>

12. TheSoundWell Vibro-therapy | vibroacoustic-training, accessed May 19, 2026, <https://www.vibroacoustic-therapy.shop/>
13. somatron-style | vibroacoustic therapy, accessed May 19, 2026, <https://www.somatron-style.com/>
14. VibroAcoustic Therapy-The Healing Power of Low Frequency Sound - healingsounds.com, accessed May 19, 2026, <https://www.healingsounds.com/vibroacoustic-therapy/>
15. TheSoundWell Products ON SALE, BUY NOW - Rehabmart.com, accessed May 19, 2026, <https://www.rehabmart.com/brand/thesoundwell.htm>
16. Benefits of Vibrational Therapy - Sound Meditation Presents, accessed May 19, 2026, <https://www.soundmeditationpresents.com/post/benefits-of-vibrational-therapy>
17. A Blend of Sound and Vibrations: What Is Vibroacoustic Therapy and How Does It Work? - Rehabmart.com, accessed May 19, 2026, <https://www.rehabmart.com/post/vibroacoustic-therapy>
18. Advances in vibration therapy for the treatment of osteoporosis - Frontiers, accessed May 19, 2026, <https://www.frontiersin.org/journals/endocrinology/articles/10.3389/fendo.2025.1611677/full>
19. Effects of broad frequency vibration on cultured osteoblasts | Request PDF - ResearchGate, accessed May 19, 2026, https://www.researchgate.net/publication/10987218_Effects_of_broad_frequency_vibration_on_cultured_osteoblasts
20. The Effects and Mechanisms of Low-Intensity Pulsed Ultrasound on Bone Remodeling: From Laboratory to Clinic - MDPI, accessed May 19, 2026, <https://www.mdpi.com/2218-273X/15/10/1351>
21. How Vibration Therapy Works | Osteoboost's Science-Driven Approach, accessed May 19, 2026, <https://www.osteoboost.com/science-of-vibration-therapy>
22. Vibration therapy and osteoporosis, accessed May 19, 2026, <https://theros.org.uk/media/uvsdna2v/vibration-therapy-fact-sheet-august-2017.pdf>
23. Possible Mechanisms for the Effects of Sound Vibration on Human Health - PMC, accessed May 19, 2026, <https://pmc.ncbi.nlm.nih.gov/articles/PMC8157227/>
24. Whole-Body Vibration Therapy for Osteoporosis: State of the Science - ACP Journals, accessed May 19, 2026, <https://www.acpjournals.org/doi/10.7326/0003-4819-155-10-201111150-00006>
25. Low-frequency vibration treatment of bone marrow stromal cells induces bone repair in vivo, accessed May 19, 2026, <https://pmc.ncbi.nlm.nih.gov/articles/PMC5243970/>
26. Whole-body vibration exercise in postmenopausal osteoporosis - PMC - NIH, accessed May 19, 2026, <https://pmc.ncbi.nlm.nih.gov/articles/PMC4440196/>
27. Osteoboost - FDA-Cleared Vibration Therapy for Low Bone Density, accessed May 19, 2026, <https://www.osteoboost.com/>
28. Osteoboost Clinical Trial | Proven to Reduce Bone Loss in Women with Osteopenia, accessed May 19, 2026, <https://www.osteoboost.com/osteoboost-clinical-trial>

29. Effects of Broad Frequency Vibration on Cultured Osteoblasts, accessed May 19, 2026,
<https://zeptoring.com/wp-content/uploads/2024/10/Studie-in-FASEB-Journal-Effect-s-of-Broad-Frequency-Vibration-on-Cultured-Osteoblasts.pdf>
30. Vibration Therapy for Osteoporosis: Does It Help or Hurt? - Healthline, accessed May 19, 2026, <https://www.healthline.com/health/vibration-plate-for-osteoporosis>
31. Clinical Benefits of Vibration Therapy [According to a Therapist] - Rehabmart.com, accessed May 19, 2026,
<https://www.rehabmart.com/post/benefits-of-vibration-therapy>
32. Advances in vibration therapy for the treatment of osteoporosis - PMC - NIH, accessed May 19, 2026, <https://pmc.ncbi.nlm.nih.gov/articles/PMC12399377/>
33. Study Details | NCT00420940 | The Influence of Vibration on Bone Mineral Density in Women Who Have Weak Bones After Menopause, accessed May 19, 2026, <https://clinicaltrials.gov/study/NCT00420940>
34. The Soundwell Body Pillow Yogibo Max for Vibroacoustic Therapy with Vibrating Pads and Full Range Music Speaker - Rehabmart.com, accessed May 19, 2026,
<https://www.rehabmart.com/product/body-pillow-thesoundwell-yogibo-max-53702.html>
35. Vibroacoustic Sound Therapy Contraindications - Sanctuary Float Center, accessed May 19, 2026,
<https://sanctuaryfloatcenter.com/vibroacoustic-sound-therapy-contraindications/>
36. Vibroacoustic Water bed Safety and Cations - Vibrating Waterbed, accessed May 19, 2026,
<https://vibroacousticwaterbed.com/vibroacoustic-water-bed-safety-and-cations/>
37. Vibroacoustic Therapy - Revive Spa, accessed May 19, 2026,
<https://reviveselfspa.com/treatment/vibroacoustic-therapy/>